StypCel™

Absorbable Hemostat

INSTRUCTIONS FOR USE

Contents

Instructions for Use Absorbable Hemostat (StypCel™)	P1
Gebrauchsanweisungen Absorbierbares Hämostat (StypCel™)	P6
Instrucciones de Uso Hemostático absorbible (StypCel™)	P12
(StypCel™) تعليمات الاستخدام مضاد النزف القابل للامتصاص	P17

Instructions for Use Absorbable Hemostat


Description



StypCel™ is a non-woven fabric which is based on absorbable material prepared by the oxidation of regenerated cellulose. This non-woven hemostat consists of multilayer of oxidized regenerated cellulose fiber, appearing white or slight yellow. It offers more convenience than the knitted fabric, because it allows the surgeon to grasp with any amount of StypCel™ that the bleeding site needs. A slight discoloration may occur over time, but it does not impact the performance. The oxidized regenerated cellulose is widely used as hemostat in clinic for decades. And it degrades or disappears within a period of about 7-14 days judged by visual inspection.

Action mechanism and degradation

Oxidized regenerated cellulose presents multiple mechanisms of action, including physical and mechanical actions in tamponade, and gel formation, and then surface interactions with proteins, platelets, intrinsic and extrinsic pathway activation. The mechanism of action whereby StypCel™ Absorbable Hemostat accelerates clotting, it is a physical effect rather than any alteration of the normal physiologic clotting mechanism. After StypCel™ Absorbable Hemostat has been saturated with blood, it swells into a gelatinous mass which aids in the formation of a clot, thereby serving as a hemostatic adjunct in the control of local hemorrhage. Absorption depends upon several factors including the amount used, degree of saturation with blood, and the tissue bed. When used properly in minimal amounts, StypCel™ Absorbable Hemostat is absorbed from the sites of implantation with practically minimal tissue reaction. It is expected to achieve hemostasis within 5min. After achieving hemostasis, the residue product can be degraded into glucose and/or glucuronic acid to be absorbed absolutely by body within 3~6 months.

Models and Specifications

Product Shape	Model
<p>Square or Rectangular</p> 	<p>HF-A-116s, HF-A-216s, HF-A-226s, HF-A-636s, HF-A-666s, HF-A-726s, HF-A-766s, HF-A-776s, HF-A-826s, HF-A-866s, HF-A-886s, HF-A-926s, HF-A-966s, HF-A-986s, HF-A-996s, HF-A-426s, HF-A-446s</p> <p>HF-A-116, HF-A-216, HF-A-226, HF-A-636, HF-A-666, HF-A-726, HF-A-766, HF-A-776, HF-A-826, HF-A-866, HF-A-886, HF-A-926, HF-A-966, HF-A-986, HF-A-996, HF-A-426, HF-A-446</p> <p>HF-A-118s, HF-A-218s, HF-A-228s, HF-A-638s, HF-A-668s, HF-A-728s, HF-A-768s, HF-A-778s, HF-A-828s, HF-A-868s, HF-A-888s, HF-A-928s, HF-A-968s, HF-A-988s, HF-A-998s, HF-A-428s, HF-A-448s</p> <p>HF-A-118, HF-A-218, HF-A-228, HF-A-638, HF-A-668, HF-A-728, HF-A-768, HF-A-778, HF-A-828, HF-A-868, HF-A-888, HF-A-928, HF-A-968, HF-A-988, HF-A-998, HF-A-428, HF-A-448</p> <p>HF-A-1110s, HF-A-2110s, HF-A-2210s, HF-A-6310s, HF-A-6610s, HF-A-7210s, HF-A-7610s, HF-A-7710s, HF-A-8210s, HF-A-8610s, HF-A-8810s, HF-A-9210s, HF-A-9610s, HF-A-9810s, HF-A-9910s, HF-A-4210s, HF-A-4410s</p> <p>HF-A-1110, HF-A-2110, HF-A-2210, HF-A-6310, HF-A-6610, HF-A-7210, HF-A-7610, HF-A-7710, HF-A-8210, HF-A-8610, HF-A-8810, HF-A-9210, HF-A-9610, HF-A-9810, HF-A-9910, HF-A-4210, HF-A-4410</p>

Product Shape	Model
<p>Round or Oval</p> 	<p>HF-B-118s, HF-B-218s, HF-B-228s, HF-B-428s HF-B-118, HF-B-218, HF-B-228, HF-B-428</p>
<p>Sector</p> 	<p>HF-C-218s, HF-C-428s, HF-C-218, HF-C-428</p>

Intended Purpose

This product is used to assist in the control of capillaries, veins, and small artery bleeding when surgical ligation and conventional control are impractical or invalid in surgery.

Indications

StypCel™ is used for assisting in the control of capillaries, veins, and small artery bleeding in traumatic hemorrhage, angiorrhagia and other lesions lead to bleeding in various surgical procedures when surgical ligation and conventional control are impractical or invalid.

Patient Target Groups

The product can be used for patients with oozing bleeds in capillaries, veins and small artery at the target area in traumatic hemorrhage, angiorrhagia and other lesions lead to bleeding in various surgical procedures.

The safety and effectiveness of StypCel™ in surgery has been studied in adult patients with no restriction of gender. The patients should not have the history of severe allergy and serious immunodeficiency and their health status need to meet the standard requirements for the surgical procedures.

Intended User& Environment

The product must be used by clinicians, and they must carefully read the instruction before use.

Operation in surgery room with professional health care facilities. Normative and local regulations must be observed.

Contraindications

- 1.StypCel™ should not be implanted in bone defects, because it may interfere with callus formation and there is a possibility of cyst formation.
- 2.StypCel™ should not be used in the control of hemorrhage from large arteries.
- 3.StypCel™ should not be used as an anti-adhesion product.
- 4.StypCel™ should not be used on the non-hemorrhagic serous oozing surface, since only the whole blood, can react with StypCel™ to produce satisfactory hemostatic effect.
- 5.Stuffing and packing is sometimes medically necessary, but the hemostat should not be used in this manner, unless it is removed after hemostasis is achieved.

Warnings

- 1.StypCel™ is a sterile medical device for single use and should not be re-sterilized or reused. Reuse of the device might pose a risk of infection to patients or hemostasis failure. The product is not suitable for autoclave or ethylene oxide sterilization.
- 2.This product acts as an auxiliary for hemostasis, it is not a substitute for suture or ligation.

3.It should be avoided to close the hemostat in a contaminated wound without drainage because it could cause complications.

4.Although the hemostat may be left in situ, it is advisable to remove the hemostat once hemostasis is achieved no matter what type of surgery. Because it has been reported in the procedures such as lobectomy, laminectomy and repair of skull fracture and lacerated lobe, that the material migrated from the site of application into foramina in bone around the spinal cord resulting in paralysis, and another case reported that the hemostat caused blindness when left into orbit of the eye. Meanwhile the material must be removed after hemostasis when used in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm regardless of the type of surgical procedure because the product, by swelling, may exert pressure resulting in paralysis or nerve damage.

5.StypCel™ should be used in dry state and should not be used with a wet instrument because it will stick to the instrument. And any wetting will weaken the hemostatic effect.

6.StypCel™ should not be impregnated with anti-infective agent or with other materials such as buffering, thrombin, because the activity of which is destroyed by the acidity of the product.

Precautions

1.According to the bleeding condition of the wound, the minimum amount of hemostatic materials should be taken out and placed on the bleeding site, and a certain pressure should be applied. In order to promote absorption and minimize the possibility of foreign body reactions, the excess hemostatic materials should be removed as far as possible before closing the surgical closure, and the removal operation should be carried out to ensure that there is no risk of re-bleeding.

2.Because the hemostat absorption could be prevented in chemically cauterized area, the StypCel™ should not be used when the wound has been treated by silver nitrate or other escharotic chemicals.

3.Make sure StypCel™ is not inhaled by patient when in the surgery.

4.If StypCel™ is temporarily placed in the cavity of a large open wound, it should not overlap the skin edges and should be removed by irrigation with sterile water or saline solution after hemostasis is achieved.

5.Surgical instruments might be used during the insertion or placement of the device, including surgical scissors, tweezers, clamp and bipolar electrocoagulation tweezers, etc.

6.If the sterile barrier system of StypCel™ is damaged, please do not use the product, and dispose it as medical waste. Because the product may be contaminated and cause failure in product performance.

7.As the storage time increases, the color of the StypCel™ may change from white to light yellow, but the performance of StypCel™ will still meet the requirements of the product.

8.StypCel™ should be used in dry state because wetness will affect the blood absorption while in use.

9.If the product appears fragmentation, badly discolored or foreign matter contaminated, please do not use the product, and dispose it as medical waste.

10.If hemostasis fails, therapeutic alternatives, such as gelatin-based hemostatic agents mixed with thrombin and fibrin sealants, or maneuver compression, can be considered.

11.Care should be taken not to apply StypCel™ too tightly when it is used as a wrap during vascular surgery.

12.Precautions should be taken in otorhinolaryngologic surgery to assure that none of the material is aspirated by the patient. (Examples: controlling hemorrhage after tonsillectomy and controlling epistaxis.)

13.StypCel™ should be used and care must be exercised to prevent plugging of the urethra, ureter, or a catheter by dislodged portions of the product.

14.Use with caution in severely infected areas that have not been effectively controlled.

15.Avoid use in patients allergic to regenerated oxidized cellulose.

Adverse reaction

1. When the hemostat was used around, in, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/ or the optic nerve and chiasm, the material left in situ may cause paralysis and nerve damage.
2. Blindness has been reported in connection with surgical repair of a lacerated left frontal lobe when the hemostat was placed in the anterior cranial fossa.
3. In the procedures such as lobectomy, laminectomy and repair of skull fracture and lacerated lobe, the material may migrate from the site of application into foramina in bone around the spinal cord resulting in paralysis.
4. In the anterior cervical arthrodesis, material could migrate from the application area to the throat resulting in respiratory obstruction.
5. Possible complications that may occur with any hemostatic used in surgery include, but not limited to, fever, infection, inflammation, hematoma, granuloma and allergy. Respiratory obstruction or pulmonary embolism, compression of the spinal cord, hearing loss, paralyzation, epilepsy may occur when used in neurosurgery.
6. Possible prolongation of drainage in cholecystectomies and difficulty passing urine per urethra after prostatectomy have been reported. There has been reported a blocked ureter after kidney resection, in which postoperative catheterization was required.

Directions for Use

1. Inspect the package and read the Instructions for Use carefully prior to use.
2. Tear the package and take out the primary package by sterile operation. The primary package is sterile and can be placed directly in a sterilized environment.
3. It is recommended that the selected size should fully cover the margins of the bleeding wound. The choice of appropriate layers depends on the surgeon's diagnosis of severity of bleeding, preferably using allowable minimal amount to achieve hemostasis.
4. Lay StypCel™ on the bleeding site, press and hold it firmly until hemostasis is achieved. Remove any excess before surgical closure in order to facilitate absorption and minimize the possibility of foreign body reaction.
5. The remaining material should be disposed as medical waste.

Tips: StypCel™ can be cut by scissors, and grasped by tweezers or forceps to the bleeding sites and pressed. And the extra hemostat can be removed by tweezers or forceps and suction system after achieving hemostat. It's allowed to perform point coagulation by electrocoagulation tweezers or forceps without removing material from the wound surface.

Expected Clinical Benefits

The clinical benefits of StypCel™ include but are not limited to:

- Effectively assist the control of capillaries and veins and small artery bleeding in traumatic hemorrhage, angiorrhagia and other lesions lead to bleeding in various surgical procedures when surgical ligation and conventional control are impractical or invalid.
- Avoiding the need of further surgical procedures by degrading completely within 3~6 months.

Material and Substances

The material and substances can be exposed to patients, including:

- Oxidized regenerated cellulose ($\leq 3\text{g/piece}$);
- 2,2,6,6-Tetramethyl-1-piperidinyloxy ($\leq 0.05\%$);
- Total nitrogen residual ($\leq 0.5\%$);
- Total chlorine and bromide residue ($\leq 1\mu\text{g/mL}$);
- Total heavy metal content ($\leq 1\mu\text{g/mL}$);

Sulfated ash ($\leq 1.5\%$);

Ethanol ($\leq 0.5\%$).

Sterilization

StypCel™ is sterilized by gamma radiation and supplied as sterile.

Shelf Life

36 months.

Storage

This product should be stored in clean and ventilated environment between 15°C and 30°C, free from corrosive gas.

Symbols used on labeling

	: Do not re-use		: Do not re-sterilize
	: Caution		: Sterilized using irradiation
	: Batch code		: Consult instructions for use or consult electronic instructions for use
	: Catalogue number		: Do not use if package is damaged and consult instructions for use
	: The storage temperature limit is 15°C to 30°C.		: Manufacturer
	: CE marking of conformity and identification number of Notified Body		: Use-by date
	: MR safe		: Single sterile barrier system with protective packaging inside
	: Medical device		: Unique device identifier
	: Date of manufacture		

The Link of Summary of Safety and Clinical Performance (SSCP)

The SSCP is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI (4260333310003PJ). The URL to the Eudamed public website is <https://ec.europa.eu/tools/eudamed>.

For the electronic version of instruction for use (e-IFU), please refer to the website: www.medprin.com.

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Any serious incident related to this device, please report to the manufacturer (Email: ae@medprin.com) and the competent authority of the Member State.

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


Form	Modell
Quadrat oder Rechteck 	HF-A-116s, HF-A-216s, HF-A-226s, HF-A-636s, HF-A-666s, HF-A-726s, HF-A-766s, HF-A-776s, HF-A-826s, HF-A-866s, HF-A-886s, HF-A-926s, HF-A-966s, HF-A-986s, HF-A-996s, HF-A-426s, HF-A-446s HF-A-116, HF-A-216, HF-A-226, HF-A-636, HF-A-666, HF-A-726, HF-A-766, HF-A-776, HF-A-826, HF-A-866, HF-A-886, HF-A-926, HF-A-966, HF-A-986, HF-A-996, HF-A-426, HF-A-446 HF-A-118s, HF-A-218s, HF-A-228s, HF-A-638s, HF-A-668s, HF-A-728s, HF-A-768s, HF-A-778s, HF-A-828s, HF-A-868s, HF-A-888s, HF-A-928s, HF-A-968s, HF-A-988s, HF-A-998s, HF-A-428s, HF-A-448s HF-A-118, HF-A-218, HF-A-228, HF-A-638, HF-A-668, HF-A-728, HF-A-768, HF-A-778, HF-A-828, HF-A-868, HF-A-888, HF-A-928, HF-A-968, HF-A-988, HF-A-998, HF-A-428, HF-A-448 HF-A-110s, HF-A-210s, HF-A-220s, HF-A-630s, HF-A-660s, HF-A-720s, HF-A-760s, HF-A-770s, HF-A-820s, HF-A-860s, HF-A-880s, HF-A-920s, HF-A-960s, HF-A-980s, HF-A-990s, HF-A-420s, HF-A-440s HF-A-110, HF-A-210, HF-A-220, HF-A-630, HF-A-660, HF-A-720, HF-A-760, HF-A-770, HF-A-820, HF-A-860, HF-A-880, HF-A-920, HF-A-960, HF-A-980, HF-A-990, HF-A-420, HF-A-440

Form	Modell
Kreis oder Oval 	HF-B-118s, HF-B-218s, HF-B-228s, HF-B-428s HF-B-118, HF-B-218, HF-B-228, HF-B-428
Kreisektor 	HF-C-218s, HF-C-428s, HF-C-218, HF-C-428

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

Forma del producto	Modèle
<p>Cuadrado o rectángulo</p> 	<p>HF-A-116s, HF-A-216s, HF-A-226s, HF-A-636s, HF-A-666s, HF-A-726s, HF-A-766s, HF-A-776s, HF-A-826s, HF-A-866s, HF-A-886s, HF-A-926s, HF-A-966s, HF-A-986s, HF-A-996s, HF-A-426s, HF-A-446s</p> <p>HF-A-116, HF-A-216, HF-A-226, HF-A-636, HF-A-666, HF-A-726, HF-A-766, HF-A-776, HF-A-826, HF-A-866, HF-A-886, HF-A-926, HF-A-966, HF-A-986, HF-A-996, HF-A-426, HF-A-446</p> <p>HF-A-118s, HF-A-218s, HF-A-228s, HF-A-638s, HF-A-668s, HF-A-728s, HF-A-768s, HF-A-778s, HF-A-828s, HF-A-868s, HF-A-888s, HF-A-928s, HF-A-968s, HF-A-988s, HF-A-998s, HF-A-428s, HF-A-448s</p> <p>HF-A-118, HF-A-218, HF-A-228, HF-A-638, HF-A-668, HF-A-728, HF-A-768, HF-A-778, HF-A-828, HF-A-868, HF-A-888, HF-A-928, HF-A-968, HF-A-988, HF-A-998, HF-A-428, HF-A-448</p> <p>HF-A-1110s, HF-A-2110s, HF-A-2210s, HF-A-6310s, HF-A-6610s, HF-A-7210s, HF-A-7610s, HF-A-7710s, HF-A-8210s, HF-A-8610s, HF-A-8810s, HF-A-9210s, HF-A-9610s, HF-A-9810s, HF-A-9910s, HF-A-4210s, HF-A-4410s</p> <p>HF-A-1110, HF-A-2110, HF-A-2210, HF-A-6310, HF-A-6610, HF-A-7210, HF-A-7610, HF-A-7710, HF-A-8210, HF-A-8610, HF-A-8810, HF-A-9210, HF-A-9610, HF-A-9810, HF-A-9910, HF-A-4210, HF-A-4410</p>

Forma del producto	Modèle
<p>Redondo u oval</p> 	<p>HF-B-118s, HF-B-218s, HF-B-228s, HF-B-428s HF-B-118, HF-B-218, HF-B-228, HF-B-428</p>
<p>Sector</p> 	<p>HF-C-218s, HF-C-428s, HF-C-218, HF-C-428</p>

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الطراز	شكل المنتج
HF-A-116s, HF-A-216s, HF-A-226s, HF-A-636s, HF-A-666s, HF-A-726s, HF-A-766s, HF-A-776s, HF-A-826s, HF-A-866s, HF-A-886s, HF-A-926s, HF-A-966s, HF-A-986s, HF-A-996s, HF-A-426s, HF-A-446s	
HF-A-116, HF-A-216, HF-A-226, HF-A-636, HF-A-666, HF-A-726, HF-A-766, HF-A-776, HF-A-826, HF-A-866, HF-A-886, HF-A-926, HF-A-966, HF-A-986, HF-A-996, HF-A-426, HF-A-446	
HF-A-118s, HF-A-218s, HF-A-228s, HF-A-638s, HF-A-668s, HF-A-728s, HF-A-768s, HF-A-778s, HF-A-828s, HF-A-868s, HF-A-888s, HF-A-928s, HF-A-968s, HF-A-988s, HF-A-998s, HF-A-428s, HF-A-448s	
HF-A-118, HF-A-218, HF-A-228, HF-A-638, HF-A-668, HF-A-728, HF-A-768, HF-A-778, HF-A-828, HF-A-868, HF-A-888, HF-A-928, HF-A-968, HF-A-988, HF-A-998, HF-A-428, HF-A-448	
HF-A-1110s, HF-A-2110s, HF-A-2210s, HF-A-6310s, HF-A-6610s, HF-A-7210s, HF-A-7610s, HF-A-7710s, HF-A-8210s, HF-A-8610s, HF-A-8810s, HF-A-9210s, HF-A-9610s, HF-A-9810s, HF-A-9910s, HF-A-4210s, HF-A-4410s	
HF-A-1110, HF-A-2110, HF-A-2210, HF-A-6310, HF-A-6610, HF-A-7210, HF-A-7610, HF-A-7710, HF-A-8210, HF-A-8610, HF-A-8810, HF-A-9210, HF-A-9610, HF-A-9810, HF-A-9910, HF-A-4210, HF-A-4410	

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